Pediatric Heart Network: Infant Single Ventricle Trial

Form S320: B-Type Natriuretic Peptide (BNP) Specimen Form

# INSTRUCTIONS: This form should be completed by site Study Coordinator and entered in Adept. Complete Form S320 for each blood sample collected for Serology Core Laboratory analysis and the file in the patient's folder.

## Section A: KEY IDENTIFYING INFORMATION

A1.	Study	Identification	Number
-----	-------	----------------	--------

#### Replaced by blinded ID

blind\_id Blinded ID

A2. Acrostic Identifier

#### Removed to protect privacy

A3. Study visit

STUDY VISIT 3 (Pre-Glenn) ......3

# STUDY VISIT 6 (Age 14 mo) ......6

 $\overline{M} = \overline{M} / \overline{D} = \overline{D} / \overline{Y} = \overline{Y} = \overline{Y} = \overline{Y}$ 

ISIT A3. Visit
----------------

A4. Date form completed

# Replaced by age at form completion

comp_age	A4. <created var="">Age at form completion,</created>
	days

A5. Name of person completing form

PRINT FULL NAME

INITIALS

### Removed to protect privacy

# Section B: SPECIMEN INFORMATION

# PATIENT SHOULD BE PLACED IN A QUIET ROOM FOR 15 MINUTES IN A SITTING OR SUPINE POSITION BEFORE THE VENOUS SAMPLE IS DRAWN.

B1. Specimen collection date

$$\frac{1}{\mathsf{M}} \cdot \frac{\mathsf{M}}{\mathsf{M}} \cdot \frac{\mathsf{M}}{\mathsf{D}} \cdot \frac{\mathsf{M}}{\mathsf{D}} \cdot \frac{\mathsf{M}}{\mathsf{M}} \cdot$$

### Replaced by age at specimen collection

bspec_age	B1. <created var="">Age at specimen collection,</created>
	days

BNP Specimen Form	Form S320	Version B: 06-05-06	Page 1 of 2
			<b>o</b>

Pediatric Heart Network: Infant Single Ventricle Trial					PHN-02
Form S320: B-Type Natriuretic Peptide (BNP) Specimen Form					111102
B2. Patient	state at blood draw	CALM	1 A	GITATED	2
PTSTATE	B2. Patient state at blood draw				

B3. Date sample shipped to Core Lab			1		1			
b3. Date sample shipped to core cab		· *	/		/			
[Not data entered]	М	Μ	D	D	Y	Y	Y	Y
			_	_	•			•

BNP Specimen Form Form S320 Version B: 06-05-0	Page 2 of 2
--	-------------